# IMPORTANT SAFEGUARDS

When using your electrical appliance, basic precautions should always be followed, including the following:

**READ ALL INSTRUCTIONS AND SAFEGUARDS BEFORE USING THIS MASSAGER.**

## DANGER

**To reduce the risk of death or injury by electric shock:**

1. Do not reach for an appliance that has fallen into water. Unplug immediately.
2. Do not use while bathing or in the shower.
3. Do not place or store appliance where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid.
4. Always unplug this appliance from the electrical outlet immediately after using and before cleaning.
5. Do not use massager without attachments.

## WARNING

**To reduce the risk of burns, electrical shock or injury to persons:**

1. This appliance should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
2. Do not operate this appliance under a blanket or pillow. Excessive heating can occur and cause fire, electric shock or injury to persons.
3. This appliance can be used by children aged from 14 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction, by a person responsible for their safety, concerning the use of the appliance in a safe way and understand the hazards involved. Children should be supervised to ensure that they do not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
4. This massager contains an electro magnetic motor. To avoid risk of injury, it must not be used within six inches of an implantable cardiac device.
5. Use this appliance only for its intended use, as described in this manual. Do not use attachments not recommended by the manufacturer.
6. Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return to Wahl for examination and repair.
7. Do not carry this appliance by supply cord or use cord as a handle.
8. Keep the cord away from heated surfaces.
9. Never drop or insert any object into any opening on this appliance.
10. Do not use outdoors.
11. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
12. Consult your physician before using this massager in case of pregnancy, varicose veins, recent surgery, unexplained pain, or illnesses including but not limited to heart condition, diabetes, or osteoporosis.
13. Do not use massager on:
  - a. Sensitive skin, swollen or inflamed areas, or skin eruptions;
  - b. Areas of poor circulation, or insensitve skin;
  - c. The front of the neck;
  - d. Genital areas;
  - e. Sleeping or unconscious persons.
14. When using a massage attachment, make sure it is properly secured onto the stem, and keep it on while massaging.
15. Do not use damaged attachments. Discontinue use of attachments if cracked or broken. May cause injury to skin.
16. Do not use this massager for more than 20 minutes in any one area.
17. To disconnect, turn all controls to the OFF position, then remove plug from outlet.
18. Do not use massager in close proximity to loose clothing or jewelry.
19. Keep long hair away from massager while in use.

**Some of our massagers have heat attachments and gel packs.**

**For a unit that has heat:**

20. Use heated surfaces carefully. May cause serious burns. Do not use over sensitive skin areas or in the presence of poor circulation. The unattended use of heat by children or incapacitated persons may be dangerous.

**For a unit that has a gel pack:**

21. Never apply the frozen gel pack directly to skin as direct application may cause frostbite. Never apply the heated gel pack directly to skin as direct application may cause burns.
22. The gel pack and contents are safe when used as directed. If the gel pack is accidentally punctured or tears, discontinue using immediately. If the gel comes in contact with the skin, wash with soap and water. If skin becomes irritated or if the gel

enters the eyes or mouth, contact a physician immediately. Keep gel pack out of reach of children and pets to avoid risks posed if gel pack is punctured or torn.

23. Check local laws and regulations for proper disposal of the gel pack.

## SAVE THESE INSTRUCTIONS FOR HOUSEHOLD USE ONLY

To reduce the risk of electric shock, this appliance has a polarized plug (one blade is wider than the other). This plug will fit in a polarized outlet only one way. If the plug does not fit fully into the outlet, reverse the plug. If it still does not fit, contact a qualified electrician to install the proper outlet. Do not change the plug in any way.

## CLEANING/STORAGE

**Always unplug massager before cleaning.**

To clean, simply wipe unit with a dry, soft cloth. Do not immerse in water or use any harsh cleaners or solvents. When storing unit, make sure cord is not twisted or kinked.

## SERVICE AND REPAIR

This unit has no user-serviceable parts. If you have problems with your massager, **do not return it to store where purchased.** Send it prepaid with a note explaining the problem to:

**Wahl Clipper Corporation  
Attention: Repair Department  
3001 Locust Street  
Sterling, Illinois 61081  
Telephone: (800) 767-9245**

For products purchased outside the United States, contact your Wahl distributor.

## 2-YEAR LIMITED WARRANTY

**VALID IN U.S.A. ONLY**

If your Wahl product fails to operate satisfactorily within two (2) years from the date of original purchase or receipt as a gift, Wahl will, at our option, repair or replace the product and return it to you FREE OF CHARGE. **Do not return this product to store where purchased.** Products needing repair under this warranty should be returned with postage pre-paid to the following address or call 1-800-767-9245 for assistance: **Wahl Clipper Corporation, Attn: Repair Department, 3001 Locust Street, Sterling, Illinois 61081.**

When mailing, please include a note describing the defect and the date of original purchase or receipt as a gift. Products out of warranty will be repaired at our standard repair charges.

This warranty does not cover products damaged by:

1. Modification or repair by anyone other than Wahl Clipper Corporation or a Wahl authorized service center;
2. Commercial (non-consumer) use or unreasonable use.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state. No other written express warranty is given. No responsibility is assumed for incidental or consequential damages of any type or for the use of unauthorized attachments. Some states may not allow the exclusion or limitation of incidental or consequential damages, so the limitation may not apply to you. Warranty valid in U.S.A. only. Outside U.S.A., check with store for Wahl authorized distributor's warranty.

## REIMBURSEMENT

This unit may qualify for reimbursement through your work FSA (flexible spending account), Medicare, Medicaid, and some insurance companies.

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