


Wahl® Hot-Cold Shiatsu Massage Support Quickstart Guide



Plug the power cord into an AC 120-V outlet.



If hot or cold is desired, heat or freeze the two gel packs. *See SPECIAL INSTRUCTIONS.*



Place gel packs in the cover/case and fasten to the front of the massage support.



To turn the unit on, press the power button once.



This unit has a dual rotation feature. The rollers rotate in the opposite direction automatically after about 45 seconds.



To turn the unit off, hold the power button for 3 seconds. This unit automatically turns off after 20 minutes.

Wahl® Hot-Cold Shiatsu Massage Support Key Features

6 DEEP-KNEADING ROTATING

NODES imitate the pressure and kneading motions of a shiatsu massage to deliver powerful relief for neck, back, shoulder, abdomen, legs, feet pain.

NODES AUTOMATICALLY


CHANGE DIRECTIONS every 45 seconds for a thorough massage.

THE ATTACHED ELASTIC

STRAP holds the cushion in place. Make sure the unit is flush with the chair or seat to allow for correct massage ergonomics.

ADD HEAT for a gentle warmth to enhance blood circulation and soothe aching muscles.

ADD COLD to aid acute and chronic pain, muscle pain, arthritis, inflammation, swelling, and injury recovery.



If you are not satisfied with your Wahl Hot-Cold Shiatsu Massage Support, please feel free to contact us at massagers.wahl.com or 1-800-767-9245.