
Special Instructions for Use of Hot-Cold Massage Therapy

DEEP HEAT ATTACHMENT

NOTE: Only the deep heat attachment can be used for the heat function.

To insert Deep Heat Attachment, simply align the special track bar with the grooved section in the massager stem (see diagram), then press firmly onto stem. To detach, use thumb and two fingers and pull away from stem.

When Deep Heat Attachment is correctly aligned and pushed on the stem, the electrical heat element is automatically connected. Test the surface for a comfortable temperature. The massager will heat to the desired level within four minutes.

NOTE: When using heat, always check temperature before placing massager against skin. Every person reacts differently to certain temperatures. We suggest the following before use — check temperature by lightly touching with fingertips the silver heat plate on the massager. If uncomfortable to your touch, reduce temperature by (1) rubbing massaging surface on soft cloth or (2) turning off unit for a short period before applying to your skin. Keep the massager moving constantly any time heat is on. Never use heat on swollen or inflamed areas, or on skin eruptions.


NEVER use the heat applicator if the plastic housing is cracked or broken, or if the metal plate is loose, bent or lifted.

CAUTION: The two holes in the stem of the massager are like an electrical plug — electrical power is transmitted through them. **NEVER** allow water near the stem, and never insert anything into the holes.

COLD GEL ATTACHMENT

****DO NOT MICROWAVE. FOR COLD THERAPY USE ONLY****

Place gel pack in the freezer. After 1-2 hours, the gel pack will reach the desired cold temperature. It is ideal to allow the gel pack to freeze overnight for optimal cold therapy. Cold will last for 20-30 minutes then remain cool. If the gel pack is uncomfortable during use, discontinue use immediately.

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